


No. 220 High Speed Air Turbine Laboratory Handpiece System



TOTALLY NEW Vibration Free and Super Quiet!

Powerful, Vibration Free, Super Quiet, and a Smooth Cutting Performance in a Dependable and Economical Package!



- **Decades of Success**

Time-proven workhorse air turbine in use in thousands of laboratories around the world

- **High-Powered Cutting Action**

350,000+ RPM provides powerful cutting action with standard FG burs

- **Ergonomic Design**

Small, light-weight aluminum design is easy to handle with less hand fatigue

- **Rugged Turbine Cartridge**

Durable friction Cartridge/Chuck system features long term reliable performance with gross trimming and coarse cutting tasks

- **Simplified Maintenance Requirements**

New Insertion and Removal Tool Sets simplify Cartridge and Metal Chuck maintenance requirements

- **Small Size**

4.25" Long x .68" Diameter x 2.2 oz.

- **Ordering Information**



No. 220 Complete System
Also Includes Filter/Regulator,
Lubricator, Tubing and Oil

<i>Item #</i>	<i>Description</i>
01000	No. 220 Air Turbine Complete System (Includes Handpiece, Regulator, Lubricator, Foot Control, Tubing and Oil)
01108	No. 220 Air Turbine Handpiece Only
01116	Replacement Turbine (Cartridge) for No. 220
01120	Replacement Metal Chuck for No. 220 (and MaxAir)
01122	Installation (Compression) Tool for Metal Chuck
01121	Standard Removal Tool for Metal Chuck
01125	Pro Removal Tool Set for Metal Chuck (Hardened Screw + Steel Punch Set)
01109	Turbine Oil, Low Viscosity, 1 Oz. Bottle
01127	Turbine Oil, Low Viscosity, 1 Qt. Bottle

Buffalo Dental Manufacturing Co., Inc.,

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA

Tel: 516.496.7200 or 800.828.0203 **Fax:** 516.496.7751

www.buffalodental.com

S01000 1109

***Please order today from
your favorite Dealer!***