

ALIGNER THERMAL PLIERS

Hot Tipped Thermal Pliers



ALIGNER Thermal Pliers will reheat the critical areas of dental aligners and tooth positioners without causing adjacent plastic to distort. Fit and form are very important when dealing with aligners, as they function as a guide to proper tooth position, so trust in a tool to get the job done right!

ALIGNER Thermal Pliers can also be used in orthodontic appliances and mouthguards.



Item #	Description
80610	ALIGNER Thermal Plier, 110 V

MADE IN USA



Contact your dealer today, or jeannieb@mailbdm.com and we will be happy to assist you!

Buffalo Dental Manufacturing Co., Inc.
159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
Email: jeannieb@mailbdm.com
www.buffalodental.com