

# Dental Operator



Buffalo Dental Manufacturing Company, Inc.  
[www.buffalodental.com](http://www.buffalodental.com)



## Amalgam Well

Stainless steel well features a deep cone-shape and internal scoring to promote easier pick-up of amalgam.

- Standard well utilizes permanently attached foam rubber base to minimize movement

Item #	Description
87400	Standard Amalgam Well



Standard Amalgam Well

## Aprons, Patient Protection Aprons

Soft, disposable, waterproof blue plastic aprons offer greater protection than paper and feature a large 26" x 30" size.

- Economical bulk packaging in rolls of 100 aprons

Item #	Description
02300	Blue Plastic Aprons, Roll of 100



Apron

## Dappen Glass

Attractive dappen glasses feature polished beveled glass construction that helps prevent material sticking.

- Economical packages of 12

Item #	Description
46500	Dappen Glass, Amber, Package of 12
46510	Dappen Glass, Blue, Package of 12
46520	Dappen Glass, Clear, Package of 12
46530	Dappen Glass, Green, Package of 12



Dappen Glass

## Glass Slabs

Polished glass mixing slabs with beveled edges promote improved mixing and safe and comfortable handling.

Item #	Description
78510	No. 2, Glass Slab (3"x4"x1/4")
78520	No. 3, Glass Slab (2 1/2"x3"x1/4")
78530	No. 4 Glass Slab (4"x6"x1/4")
78570	No. 10 Glass Slab (3"x6"x1/2")
78580	No. 10A Glass Slab (3"x6"x3/4")



No. 10 Glass Slab

## Headrest Covers

Durable headrest covers resist moisture, protect headrests, and improve patient comfort.

- Regular Size effective for protecting double-cup headrests (10 1/2" L x 9 1/2" H)
- Regular Size with Gussets feature extra fullness and are ideal for double-cup headrests (12 1/2" x 8 1/2")

- Large Size with Gussets excellent for protecting contour chair and pillow-type headrests (15 1/4" x 8 1/2")

Item #	Description
31750	Regular, Blue, Roll of 500
31790	Regular, Clear, Roll of 500
31740	Regular, White, Roll of 500
31820	Regular w/ Gussets, Clear, Box of 500
31830	Regular w/ Gussets, White, Box of 500
31850	Large w/ Gussets, Clear, Box of 400
31860	Large w/ Gussets, White, Box of 400



Headrest Covers

## Mirror, Patient Mirror

- High quality, 3 1/2" diameter mirror with distortion-free glass and swivel head promotes improved viewing and convenient use.

Item #	Description
61755	Patient Mirror, 3 1/2" Diameter



Mirror

## Napkin Holders

Durable napkin holders with strong metal chains and alligator clips for secure holding and long life.

- Nickel-plated version fully autoclavable for improved infection control
- Colored sleeve versions (14" total length) and Nickel-plated (10 1/4" total length) available

Item #	Description
51451	Nickel-plated, Autoclavable, Pkg. of 12
51591	Blue Sleeve, Pkg. of 12
51611	Red Sleeve, Pkg. of 12
51621	White Sleeve, Pkg. of 12
51451X72	Nickel-plated, Autoclavable, Pkg. of 72
51591X72	Blue Sleeve, Pkg. of 72
51611X72	Red Sleeve, Pkg. of 72
51621X72	White Sleeve, Pkg. of 72
51623X72	Yellow Sleeve, Pkg. of 72



Nickel-plated

Color-Coded

### Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Drive • Syosset, NY 11791-0678 USA  
 Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax(orders): (516) 213-6610  
 www.buffalodental.com



### Trays, Dual Arch Fluoride Trays

Soft foam, I-beam style trays conform comfortably to most arches, treat both arches simultaneously, and require less fluoride than other trays.

- Versatile 5" long trays can be easily trimmed to desired length

**Item # Description**  
56553 Fluoride Trays, Pkg. of 100



Dual Arch Fluoride Trays

### X-Ray Rapid Developer and Instant Fixative

Rapid developer and fixer designed to allow immediate development of diagnostic x-rays with excellent visual quality.

- Compatible with Kodak brand Ultra-Speed and other films
- Use for prompt and accurate assessment of endodontic, toothache, tooth fracture and other indications

- Rapid Developer Pkg. includes 110 cc liquid, powder vial, and developing jar – simply mix all liquid and powder together in jar to get up to 30 days or more of developing power
- Instant Fixative Pkg. contains 110 liquid and fixative jar- enough to fix films for up to 30 days or more
- Simplified procedure: Immerse film in Developer for 15 sec.; Rinse; Immerse in Fixative liquid for 30 sec.

**Item # Description**  
33800 Rapid Developer  
42300 Instant Fixative



X-Ray Rapid Developer and Instant Fixative

### MOLLOPLAST-B® Lab Processed Silicone Soft Reline Material- Now Microwaveable!

MOLLOPLAST-B is an industry-leading heat-processed silicone material that stays soft permanently and features decades of documented success in denture reline cases. The material bonds firmly to new or existing acrylic dentures, and stays elastic and bacterial free for years.

- New Microwaveable consistency offers updated, fast curing and added plaque inhibitors
- Convenient lab procedure for use with existing dentures or during construction of new dentures

- Indicated for improving comfort and fit of acrylic dentures, relieving pressure spots, steadying dentures, and attachment to existing dentition
- Primo® Adhesive excellent for preparing existing acrylic for permanent attachment to MOLLOPLAST-B
- Lustrol® Polish in 1:1 mixing ratio creates protective glossy finish on final MOLLOPLAST-B surfaces

**Item # Description**  
62300 MOLLOPLAST-B, 45 gm Reg. Size  
62310 MOLLOPLAST-B, 170 gm Lab Size  
58600 Lustrol, 6 ml Base and 6 ml Catalyst New 1:1 Mixing Ratio!  
73000 Primo, 15 ml Bottle and Brush



Molloplast-B

### SoftGrind™ Grinding Sleeves

The ideal abrasive system for adjusting and finishing silicone soft reline materials.

- Unique reusable HP Mandrel features a rubber head to securely hold sleeves with a friction fit
- Sharp abrasive sleeves attach easily and remove silicone quickly with superior control

- Small sleeve size (5/16" Diameter x 1/2" Length w/ 60 degree tip angle) provides precise control

**Item # Description**  
62400 SoftGrind Intro Kit, Pkg. of 10 Sleeves and 1 HP Mandrel  
62410 SoftGrind HP Mandrel, Pkg. of 1  
62420 SoftGrind Grinding Sleeves, Pkg. of 10



SoftGrind Grinding Sleeves

### Pre-Polishers

Unique pre-polishing wheels create a smooth finish on silicone relines prior to the final polishing step.

- Convenient packaging of 10 disposable Wheels and one reusable HP Mandrel

- 3M Scotch-Brite® - style texture provides excellent finishing action
- Small size (1 1/8" Diameter x 1/4" Thick) allows improved access

**Item # Description**  
62200 Pre-Polishers, Pkg. of 10 Wheels and 1 HP Mandrel



### Mollo-Cutter™ Burs

Ultra-sharp steel burs with cross-cut flutes designed for fast and smooth trimming of silicone soft reline materials and acrylic denture base.

- Resist clogging and heat generation
- Three assorted HP shapes per package

**Item # Description**  
62280 Mollo-Cutter Burs, HP Sized, Set of 3



Mollo-Cutter Burs

### Mollobrasive™ Points

Abrasive points with special, non-clogging grit designed for excellent smoothing and finishing of silicone soft reline materials.

- Three HP points per package

**Item # Description**  
62290 Mollobrasive Points, HP Sized, Set of 3



Mollobrasive Point

**Buffalo Dental Manufacturing Co., Inc.**

159 Lafayette Drive • Syosset, NY 11791-0678 USA  
Tel: (800) 828-0203 • (516) 496-7200 • Fax: (516) 496-7751 • Fax/orders: (516) 213-6610  
[www.buffalodental.com](http://www.buffalodental.com)