

Dust Containment Unit

The solution to dust and debris from handpiece grinding.



Enhance safety for you and your assistant's with the new Buffalo Dental Manufacturing **Dust Containment Unit**. The **D.C.U** for short, is an enclosed unit featuring an internal light, suction fan and filter system that will allow you to use the supplied handpiece in a more sanitary and safe environment.

Ideal for a dentist to use in a mini-lab or backroom.

Features:

- Small Size (11" wide x 12" deep x 8 1/2" high) fits easily on countertops.
- High Powered, reliable, 35,000 RPM V35 Electric Handpiece System ideal for most cutting and grinding needs.
- Hinged top window and internal light fixture provide bright illumination for clear viewing.
- Dual, wide hand holes allow comfortable hand positioning and access to large workspaces.
- Efficient suction fan and filter collect damaging grinding debris and maintain quiet operation.
- High capacity internal filter that is fully replaceable.
- Additional 120 V AC power outlet on back of the unit for operating electric devices.
- Hinged rear door for easy cleaning.
- Convenient handpiece clip on top of unit to keep your handpiece at the ready.



The Buffalo Dust Containment Unit
which comes with a V35 Electric Handpiece System



Buffalo # **Description**

36800 *Dust Containment Unit, 120 V AC*

Buffalo Dental Manufacturing Co., Inc.
159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
Email: jeannieb@mailbdm.com
www.buffalodental.com