

# ThermaKnife™

## Hot Tipped Tray Trimming Knife



***Effective and Time-Saving Knife  
For Fast and Easy Trimming of  
Bleaching Trays, Night Guards,  
Mouthguards, and Other  
Thermoplastic Materials***

ThermaKnife is the improved way to quickly and easily trim bleaching trays, mouthguards, and many other thermoplastic materials.

Simply choose the pen with the tip design best for your needs, set the console to the desired temperature, and complete your trimming job with ease. The ThermaKnife has the ability to support two

different pens at one time. By using the A and B switch on the front of the unit the user can switch between sending heat from one pen to the other!

The adjustable temperature control makes it easy to customize the process to your cutting speed

and material thickness needs, and the tip heats up and is ready to use in seconds.

Optional pens are available featuring different tip designs, the pens can be interchanged quickly and easily, and they can be stored safely in the pen rest when not in use.



### Features/Benefits

- Fully adjustable temperature control for precise cutting action
- Cutting tip heats ultra-fast to save time and eliminate waiting
- Lightweight pen features an ergonomic design and sharp, permanent tip for improved comfort and control
- A & B selection switch allows for the use of two different Thermal Pens to reduce time and increase trimming productivity
- Fast and easy pen replacement (via snap connector at back of pen) for added convenience
- Convenient pen rest provides added safety when unit not in use
- Small portable size: Console 5.75" W x 5" D x 4"H (incl. Rest)  
Pen: 6" L x .75" W (at foam grip,) 2.5 lb. total console/pen weight



**Knife-Edge Tip**



**Concave Tip**



**Waxing Tip**

### Ordering Information

Description	Item No.
ThermaKnife Intro Kit, 120V AC (Includes Knife Edge Pen)	80500
Replacement Pen with Concave Tip	80510
Replacement Pen with Knife Edge Tip	80520
Replacement Pen with Waxing Tip	80530

**Buffalo Dental Manufacturing Co., Inc.,**  
159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA  
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751  
[www.buffalodental.com](http://www.buffalodental.com)

S80500 0706

*Please order today from  
your favorite Dealer!*