

UV-MT Model Trimmer



Ultra-Violet Model Trimmer Kills Germs While It Turns!

Dental stone models contain various bacteria and organisms which can be spread throughout the lab.

Made in USA!

Powerful Internal UV-LED Diode kills Bacteria and other organisms on wheel, in grindings and on model. Reduces cross contamination from model to model.

- **Powerful and Fast Trimming**

Industry-leading extra-high torque motors and aggressive-cutting coarse Carborundum wheels cut fast and smooth saving time and money

- **Solid Construction**

Classic, heavy construction and time-proven balanced direct-drive systems run quietly and can stay steady and trouble-free on the job for decades

- **Enlarged Window**

Larger work window provides increased working area, improved vision, and full wheel use while resisting splashing

- **Versatile Table**

Calibrated work table indicates model angles and can be easily tilted and repositioned to suit nearly every trimming requirement

- **Complete Ready-to-Run System**

All systems ship complete with a pre-mounted balanced Aggressive-cutting Carborundum Trimming Wheel, ¼" ID Water Supply Tubing, 1" ID Water Exhaust Hose, and Wheel Wrenches

- **Optional Accessories**

A wide range of optional accessories are available including pre-mounted water solenoid valves, rubber feet, safety shields, numerous diamond and conventional abrasive wheels, disposable disks, and several disposable plaster trap options

- **Ordering Information**

<u>Item Number</u>	<u>Description</u>
--------------------	--------------------

61780	10" Diameter 1/3 HP Ultra-Violet Model Trimmer, 120 V AC, 60 Hz
61785	12" Diameter 1/3 HP Ultra-Violet Model Trimmer, 120 V AC, 60 Hz



***US and Foreign Patents Pending.**

Buffalo Dental Manufacturing Co., Inc.,

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA

Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751

www.buffalodental.com

S61780 1108

Please order today from your favorite Dealer!