

Model Arch Trimmer



Premium High-Torque Trimmer Provides Fast and Quiet Cutting for the Most Demanding Model Arch Dry Trimming Applications!



- **Fast and Easy Dry Trimming**

High-torque, 5400 rpm, 1/8 HP heavy-duty motor promotes quiet, fast, and nearly-effortless cutting

- **Improved Dust Control**

Enlarged 7" x 7" work table with built-in dust control ducting and 1" OD exhaust port on back maintain quiet and effective dust control- ready for attachment to your existing dust collector

- **Rapid Tool Changes**

Quick-release work table and 1/4" collet allow convenient and fast changing of 1" maximum OD cutting tools

- **Improved Safety**

Strong 4" x 6" clear polycarbonate safety shield helps protect the operator from debris, and the interlocking on-off switch automatically resets to off in case of power interruption

- **Buffalo Tough Cutting Instruments**

System ships complete with a fast cutting Buffalo **Tough** 1/2" diamond-cut tungsten carbide bur and a 3/4" coarse-grit arbor band set

- **Compact Size**

8 in. (20.4 cm) Width x 10 in. (25.4 cm) Depth x 8-1/8 in. (20.6 cm) Height x 14 lb. (6.4 kg) Weight

- **Complete Range of Products**

Please see our complete range of model trimming equipment and supplies including Wet and Dry Model Trimmers, Trimming Wheels and Belts, Knives and Burs, Safety Shields, and Disposable Plaster Traps.



- **Ordering Information**

<u>Item #</u>	<u>Description</u>
64000	Model Arch Trimmer Complete System, 120V AC, 1.2 A, including Safety Shield, Wrench Set, 1/2" Diameter Diamond-Cut Tungsten Carbide Bur (#19400), 3/4" Arbor Band Mandrel (#02925), and Coarse Grit Zirconia Arbor Bands, Pkg. of 10 (#02700)

Made in USA!

Notes: For laboratory use only. Use extreme caution and wear personal protective eyewear, gloves and equipment when operating a Model Arch Trimmer. Avoid personal contact with the spinning bur, and always unplug the Trimmer when changing burs or cleaning the housing.

Buffalo Dental Manufacturing Co., Inc.,
129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
www.buffalodental.com

S64000 0708

**Please order today from
your favorite Dealer!**