

# Heavy Duty Bench Engines



**Some Things Never Go Out of Style!  
Classic Buffalo Dental Bench Engines  
Deliver Maximum Performance in the Most  
Demanding Nail, Jewelry, Industrial,  
Laboratory and Hobbyist Applications**

**Improved**

- Work-horse strong and durable 1/5 HP, 0-20,000 rpm motors- *Made in USA!*
- All ball-bearing construction never needs lubrication
- Sturdy wide base for improved stability
- Every system is wear-tested to ensure smooth, true-running operation at cool, user-friendly temperatures
- Beautiful polished chrome housings
- Handpieces feature sealed ball bearings for cool & true running and versatile 0 – 5/32" adjustable chucks
- Variable Speed Electronic Foot Controls feature durable solid-state circuitry for long term reliable performance
- Systems include engine, adjustable arm, drive belt, 3-wire grounded plug, 14" cord, and Variable Speed Foot Control
- One year Warranty



**Highest Quality and Most  
Reliable 1/5 HP Motors!**

## • Ordering Information

Item #	Description
36919HP	No. 16 Chrome Bench Engine Complete System w/Triple-Section Regular Arm, <u>Labtex Handpiece</u> and Variable Speed Foot Control, 120V AC
36920HP	No. 16 Chrome Bench Engine Complete System w/Bantam (Short) Arm, <u>Labtex Handpiece</u> and Variable Speed Foot Control, 120V AC
36919HP101	No. 16 Chrome Bench Engine Complete System w/Triple-Section Regular Arm, <u>No. 101 Handpiece</u> and Variable Speed Foot Control, 120V AC
36920HP101	No. 16 Chrome Bench Engine Complete System w/Bantam (Short) Arm, <u>No. 101 Handpiece</u> and Variable Speed Foot Control, 120V AC
36937HP	No. 16 Chrome Bench Engine Complete System w/Triple-Section Regular Arm, <u>Labtex Handpiece</u> and Variable Speed Foot Control, <u>220V AC</u>
36939HP	No. 16 Chrome Bench Engine Complete System w/Bantam (Short) Arm, <u>Labtex Handpiece</u> and Variable Speed Foot Control, <u>220V AC</u>
36937HP101	No. 16 Chrome Bench Engine Complete System w/Triple-Section Regular Arm, <u>No. 101 Handpiece</u> and Variable Speed Foot Control, <u>220V AC</u>
36939HP101	No. 16 Chrome Bench Engine Complete System w/Bantam (Short) Arm, <u>No. 101 Handpiece</u> and Variable Speed Foot Control, <u>220V AC</u>
49300	No. 101 Heavy Duty Lab Handpiece, 0 – 45,000 RPM
49445	Labtex Lab Handpiece, 0 – 45,000 RPM

**New Systems  
w/Handpieces!**

*Made in USA!*


**Performance testing reveals cooler -  
running, quieter, higher-speed and more  
powerful than competing Engines!**



**Labtex Handpiece**



**No. 101 Handpiece**



**Complete System**

**Buffalo Dental Manufacturing Co., Inc.,**  
129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA  
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751  
[www.buffalodental.com](http://www.buffalodental.com)

S36919 0706

**Please order today from  
your favorite Dealer!**