

Buffalo Tough Carbide Lab Burs



Made-in-USA Premium Carbide Burs Provide Outstanding Cutting Speed, Durability and Value!

• Ultra-Fast Cutting

Automated manufacturing and highest-quality materials assure ultra-fast cutting and long lasting edges to save your time

• Excellent Value

USA manufacturing location and automated machining centers assure premium quality, cut out middle men, and save you money

• Optional TNT Titanium Nitride-Coated Heads

Available with premium USA-applied thick titanium nitride-coated heads for longest-lasting edges and superior cutting

Made in USA!



• Dual Cut & Diamond Cut Heads Available

Traditional dual cut heads improve performance and minimize clogging, and diamond cut versions help pulverize materials and excel at acrylic work

• Optional Dust Shields on All HP Burs

Patented dust shields available on HP burs protect handpiece bearings from damaging dust (no extra charge- add -SH to end of HP Item No.)




Popular HP Shapes



3/32" Shank	Shape:	88A Barrel	82T Cone	83E Cylinder	52C Egg/Oval	81A/61A Flame	73C Inverted Cone	53A Pear	52D Round	84T Taper	63B Flame	51A Small Flame
Dual Cut Regular #:		18000	18010	18020	18030	18050	18060	18070	18080	18100	18120	18130
Dual Cut Regular w/ TNT #:		88A-GOLD	82T-GOLD	83E-GOLD	52C-GOLD	81A-GOLD	73C-GOLD	53A-GOLD	52D-GOLD	84T-GOLD	63B-GOLD	51A-GOLD
Diamond Cut Coarse #:		18000CDM	18010CDM	18020CDM	18030CDM	18050CDM	18060CDM	18070CDM	18080CDM	18100CDM	18120CDM	18130CDM
Diamond Cut Fine #:		18000FDM	18010FDM	NA	18030FDM	18050FDM	NA	18070FDM	18080FDM	18100FDM	18120FDM	NA

Popular 3/8" Lathe Burs



1/4" Shank	Shape:	A3/8 Taper	B3/8 Cone	C3/8 Flame	D3/8 Egg	E3/8 Cylinder	F3/8 Barrel	G3/8 Ball
Dual Cut Regular #:		19130	19100	19120	19110	19150	19140	19160
Dual Cut Regular w/ TNT #:		A3/8-GOLD	B3/8-GOLD	C3/8-GOLD	D3/8-GOLD	E3/8-GOLD	F3/8-GOLD	NA
Diamond Cut Coarse #:		19130CDM	19100CDM	19120CDM	19110CDM	19150CDM	19140CDM	NA

Please check with Buffalo Dental for additional shapes and sizes.

A full range including smaller-headed HP shapes, fine-cut versions, and popular 1/8" and 1/4" head sizes (mounted on 1/8" shanks are available).

Popular 1/2" Lathe Burs



1/4" Shank	Shape:	A1/2 Taper	B1/2 Cone	C1/2 Flame	D1/2 Egg	E1/2 Cylinder	F1/2 Barrel	G1/2 Ball
Dual Cut Regular #:		19340	19300	19330	19320	19310	19350	19360
Dual Cut Regular w/ TNT #:		A1/2-GOLD	B1/2-GOLD	C1/2-GOLD	D1/2-GOLD	E1/2-GOLD	F1/2-GOLD	NA
Diamond Cut Coarse #:		19340CDM	19300CDM	19330CDM	19320CDM	19310CDM	19350CDM	NA

Buffalo Dental Manufacturing Co., Inc.,
 129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA
 Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
www.buffalodental.com

S18100 0607

Please order today from your favorite Dealer!