

DryVac™ Dust Collector

The strong dependable suction system for situations requiring high volume dust extraction. Perfect for dry model trimming. Powerful ½ HP intermittent-duty motor provides strong, high pressure 100 CFM suction. Includes motor silencer assembly and filter paper. 3" O.D. inlet. Complete with 4 Ft. vacuum hose. 16" H x 12" D, Wt. 28 lbs.

<u>Item #</u>	<u>Description</u>	<u>Item #</u>	<u>Description</u>
34000	DryVac Complete 120 V AC	BAG3	Replacement Filter Bag for DryVac
34050	DryVac Complete 220 V AC		



DryVac



DustClear PreFilter

DustClear™ PreFilter Dust Collection System

DustClear is an efficient dust collection cabinet that features innovative AirReverse™ filter cleaning technology. The non-powered PreFilter is designed for use directly behind extremely dusty operations such as dry model trimming, but ahead of existing central dust collection systems. When used in this manner as an economical pre-filtering unit, DustClear extracts nearly 100% of the waste from dusty procedures and protects down-stream dust collection equipment. Requires 80 psi compressed air connection to operate AirReverse system. Features 2 7/8" O.D. Inlet and Outlet. Optional top-mounted 1/2 HP vacuum motor converts DustClear to a high-powered, 100 CFM, stand-alone dust collector ideal for intermittent-duty use. 11" W x 14" D x 18 1/2" H x 22 lb. (w/o Motor)

<u>Item #</u>	<u>Description</u>	<u>Item #</u>	<u>Description</u>
37005	DustClear PreFilter Unit	37000	Fully Powered DustClear, 120 V AC
37008	Vacuum Motor, 120 V AC	F114	Replacement Filter Cartridge for DustClear



DustStation I

DustStation I™

DustStation I is a high performance, commercial grade dust collection system that features powerful suction, quiet operation and a unique and effective easy cleaning filtration system. 21" H x 15" W x 15" D, Wt. 73 lbs

- Powerful 1/2 HP Continuous-Duty Motor, 400 CFM Suction, Dual 3" O.D. Inlets
- Quiet Air-tight seals and quiet motor keep noise below 72 decibels
- Unique Filter Assembly with 14 Cotton Covered Mesh filters entrap particles above .5 micron
- External Shaker Handle allows fast, effective cleaning of filters

<u>Item #</u>	<u>Description</u>	<u>Item #</u>	<u>Description</u>
35050	DustStation I, 120 V AC	35052AO	Set of 12 Filter Bags for Aluminum Oxide
35060	DustStation I, 220 V AC	35052	Set of 14 Replacement Filters
35050S	Silencer for DustStation I	35040-08	Vacuum Hose (3" I.D.), 8 Ft. Section

Mobile LabStation™ with Suction

The Mobile LabStation combines the portability and utility of a mobile cart with the performance and functionality of our Workbench Suction Drawer. The unit is perfect for the small dental office or lab desiring a high-performance portable workstation with built-in suction and enough storage space for multiple lab items. Pre-wired with 6-outlet power strip and ready to plug in! Attractive and durable white laminant finish. 32"H x 20 1/2" W x 23"D, Wt. 90 lbs.

<u>Item #</u>	<u>Description</u>
41250	Mobile LabStation w/ Suction 120 V AC



Mobile LabStation

Workbench Suction Drawer

The easy to install, fully functional, small dust collection system that conveniently mounts under your work surface and extends when needed, but tucks away when not in use. The Workbench Suction Drawer includes a continuous-duty transflow blower that creates 125 CFM suction and collects fine dust particles but remains quiet during use. It also features a unique water trap for collecting precious metal grindings, and a flip-up shield for extra eye protection. 7" H x 18 1/2" W x 19" D, Wt. 23 lbs.

<u>Item #</u>	<u>Description</u>
36000	Workbench Suction Drawer w/ Suction 120 V AC
36100	Workbench Suction Drawer w/o Suction (w/ 3" O.D. port on back)
42113	Replacement Filters, Pkg. of 3



Workbench Suction Drawer

Buffalo Dental Manufacturing Co., Inc.,

129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA

Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751

www.buffalodental.com

S34000 0607

**Please order today from
your favorite Dealer!**