

M35 Compact Electric Lab Handpiece System



New High Torque Variable Speed System Features a Compact and Versatile Design

- **High Torque Cutting Power from 0 – 35,000 RPM**

New high torque handpiece and control system provide excellent cutting power and smooth control in a compact system

- **Improved Versatility**

System is designed for 120V AC operation but can be easily switched to 220V AC via setting on the underside of the Console

- **Improved Console and Controls**

Improved controls provide fully variable speed via console or foot control, simplify operation and protect the controls and handpiece

- **Excellent High and Ultra-Low Speed Operation**

System excels at providing a full range of usable operating speeds from ultra-low to 35,000 RPM

- **Heavy Duty Construction**

Components and system designed to provide reliable and trouble-free operation. One year warranty

- **Sizes**

Handpiece: 6.25 in. (159 mm) Long x 1.1 in. (28.5 mm) Back Diameter, .79 in. (20 mm) Front Diameter x 8 oz. (.23 kg) Weight

Console: 5 in. (128 mm) Width x 6.5 in. (166 mm) Depth x 3.25 in. (82 mm) Height x 2.9 lb. (1.3 kg) Weight



The Motor and Armature feature precision dynamic balancing for improved performance and less noise and vibration

- **Ordering Information**

<i>Item #</i>	<i>Description</i>
38120	M35 Electric Lab Handpiece System Complete, 120V AC (& 220V AC), including Control Console, Electric Handpiece w/HP Collet, Foot Control, Handpiece Holder and Collet Wrenches
38140	M35 Electric Lab Handpiece Only w/HP Collet Chuck
38025BS	Square Carbon Brushes for M35 Electric Handpiece, Set of 2
37805-107-3	Replacement HP Collet Chuck for M35 Handpiece, 3/32" Size
37805-COLLET3MM	Optional 3mm Collet for M35 Handpiece
37805-COLLET18	Optional 1/8" Collet for M35 Handpiece

Notes: For laboratory use only. Please use extreme caution and wear personal protective eyewear and equipment when operating the M35 system.

Distributed by:

Buffalo Dental Manufacturing Co., Inc.,
129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA
Tel: 516.496.7200 or 800.828.0203 **Fax:** 516.496.7751
www.buffalodental.com

S38120 0706

***Please order today from
your favorite Dealer!***