

Sta-Vac™ II

Versatile Vacuum Forming System



Improved Classic High Performance Vacuum Forming System that Allows Both Traditional Flat Table and Accurate Bead Basket Adaptations



Features/Benefits

- **Versatile New ConverTable™ Combination Vacuum Table and Bead Basket**
Unique new table design can be used in quick and easy flat table mode or flipped over and used in the versatile bead basket mode for infinite model positioning and improved vacuum flow
- **High Performance Electric Vacuum Forming System**
Dual-stage electric vacuum motor provides powerful suction and highly accurate adaptations- no compressed air connections required!
- **Accurate Frame and Alignment System**
Spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table
- **Heavy Duty Construction**
Solid construction (7" W x 16" D x 11" H, 32 lbs.) for maximum durability and long life
- **Convenient and Safe Design**
Unique housing design prevents motor contamination, shields heater, and promotes ergonomic operation

Ordering Information

<u>Item Number</u>	<u>Description</u>
80193	Sta-Vac II Vacuum Forming System, 120 V AC, w/ Two ConverTables
80194	Sta-Vac II Vacuum Forming System, 220 V AC, w/ Two ConverTables
80198	Replacement ConverTable for Sta-Vac II
80188	Replacement Steel Beads for Sta-Vac II (and Accu-Vac), 3 Lb. Pkg.

Note: Systems includes Sta-Vac II unit and sample range of Vacuum Forming Materials

Buffalo Dental Manufacturing Co., Inc.,
129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
www.buffalodental.com

S80193 0706

**Please order today from
your favorite Dealer!**