

Buffalo Dental Waxing Unit



Compact and Economical 3-Compartment Wax Heating Unit with Adjustable Temperature Control

Features/Benefits

- On/Off switch located on power cord to prevent wax seepage into switch mechanism
- Controlled temperatures from 110° to 140° F for convenient operation (adjustment screw on back left)
- Provides three 11/16" deep wax compartments (one is 2-1/2" x 2-7/8", and two are 1-5/8" x 1-3/8")
- Attractive design with wood sides and helpful heat indicator light on front panel



Ordering Information

<u>Item No.</u>	<u>Description</u>
86000	Waxing Unit, 120V AC
86010	Waxing Unit, 220V AC

Buffalo Dental Manufacturing Co., Inc.,
129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA
Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
www.buffalodental.com

S86000 0706

***Please order today from
your favorite Dealer!***