

Wet Model Trimmers



Classic Wet Model Trimmers Feature Outstanding Performance and Reliability

Made in USA!

- **Powerful and Fast Trimming**

Industry-leading extra-high torque motors and aggressive-cutting coarse Carborundum wheels cut fast and smooth saving time and money

- **Solid Construction**

Classic, heavy construction and time-proven balanced direct-drive systems run quietly and can stay steady and trouble-free on the job for decades

- **Enlarged Windows**

Larger work windows provide increased working area, improved vision, and full wheel use while resisting splashing

- **Versatile Tables**

Calibrated work tables indicate model angles and can be easily tilted and repositioned to suit nearly every trimming requirement

- **Complete Ready-to-Run Systems**

All systems ship complete with a pre-mounted balanced Aggressive-cutting Carborundum Trimming Wheel, ¼" ID Water Supply Tubing, 1" ID Water Exhaust Hose, and Wheel Wrenches

- **Optional Accessories**

A wide range of optional accessories are available including pre-mounted water solenoid valves, rubber feet, safety shields, numerous diamond and conventional abrasive wheels, disposable disks, and several disposable plaster trap options



Item Number	Description
-------------	-------------

120V AC Options

61790A	10" Diameter 1/3 HP Wet Model Trimmer, 120V AC, 60 Hz
61792	10" Diameter 1/2 HP Wet Model Trimmer, 120V AC, 60 Hz
61825	12" Diameter 1/3 HP Wet Model Trimmer, 120V AC, 60 Hz
61820	12" Diameter 1/2 HP Wet Model Trimmer, 120V AC, 60 Hz

220V AC Options

61810	10" Diameter 1/3 HP Wet Model Trimmer, 220V AC, 60 Hz
61792A	10" Diameter 1/2 HP Wet Model Trimmer, 220V AC, **50 Hz**
61830	12" Diameter 1/2 HP Wet Model Trimmer, 220V AC, 60 Hz

Trimming Wheels

61940	10" Diamond Trimming Wheel w/Aluminum Backing & 1" Center Hole
61950	12" Diamond Trimming Wheel w/Aluminum Backing & 1" Center Hole
61990	10" Carborundum Trimming Wheel, Coarse Grit w/1" Center Hole
62030	12" Carborundum Trimming Wheel, Coarse Grit w/1" Center Hole
61965	10" Quick-Stick™ "Wrench-Free" Disposable Silicon Carbide Abrasive Disks, Coarse Grit (w/2" Center Hole), Pkg. of 6- <i>Requires Plastic Backing Wheel + Washer</i>

Accessories

78900	Custom Clamp-On Safety Shield w/Goose-Neck (Clamps to Trimmer Base)
78400	Standard Safety Shield w/Goose (Bolts to Counter Top)
61798	Rubber Noise-Dampening Feet for Model Trimmers, Set of 4

Buffalo Dental Manufacturing Co., Inc.,

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA

Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751

www.buffalodental.com

S61790A 0707

*Please order today from
your favorite Dealer!*