

X40 Premium Electric Lab Handpiece System



New 40,000 RPM High-Torque System Updates Our Classic X35 and Provides Additional Speed and Torque



- New 40,000 RPM Handpiece**
 40k blue electric handpiece features an improved heavy-duty spindle, upgraded bearings, and more durable square carbon brushes, providing higher speeds, smoother operation, and increased reliability
- Improved Console *Made in USA!***
 The upgraded control console features dual power outlets allowing fast and easy switching between two separate handpieces as desired (one handpiece included with system), saving time-wasting tool changes. Higher-voltage output and additional overload control help protect the handpiece from damage during extended heavy grinding



- Superior Control**
 The system provides exceptional control at ultra-low rotational speeds, high torque cutting power, and fully variable speed from 0-40,000 rpm. The X40 comes complete with a fully variable speed foot control for added versatility
- Heavy Duty Construction**
 Components and system designed to provide reliable and trouble-free operation. One year warranty
- Sizes**
Handpiece: 6.25 in. (159 mm) Long x 1 1/8 in. (29 mm) Major Diameter x 5/8 in. (16 mm) Front Diameter x 8.6 oz. (244 gm) Weight
Console: 7.25 in. (184 mm) Width x 7 in. (177 mm) Depth x 3.5 in. (89 mm) Height x 6.3 lb. (2.9 kg) Weight

• Ordering Information

<i>Item #</i>	<i>Description</i>
38085	X40 Electric Handpiece System Complete, 120V AC, including Control Console, 40k Blue Handpiece w/HP Collet, Variable Speed Foot Control, Handpiece Holder, Spare Square Carbon Brushes and Collet Wrenches
38087	X40 Electric Handpiece System Complete, 220V AC, including Control Console, 40k Blue Handpiece w/HP Collet, Variable Speed Foot Control, Handpiece Holder, Spare Square Carbon Brushes and Collet Wrenches
38025CN	Replacement HP Collet Chuck for X40 Handpiece, 3/32" Size
38025BS	Replacement Square Carbon Brushes for X40, Set of 2

For laboratory use only. Please use extreme caution and wear personal protective eyewear and equipment when operating this system. Follow our manufacturer's recommendations for maximum safe operational speeds for burs

Additional chuck sizes are available for X40 Handpieces - please call for details

Buffalo Dental Manufacturing Co., Inc.,
 129 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 USA
Tel: 516.496.7200 or 800.828.0203 **Fax:** 516.496.7751
www.buffalodental.com

S38085 0706

Please order today from your favorite Dealer!