

# Dust Containment Unit <sup>\*Patented</sup>

The solution to dust and debris from handpiece grinding.



## PROTECTION FOR YOU AND YOUR STAFF!

Enhance safety for you and your assistant with the new Buffalo Dental Manufacturing **Dust Containment Unit**. The **D.C.U** for short, is an enclosed unit featuring an internal light, suction fan and filter system that will allow you to use the **included V35 Handpiece System** in a more sanitary and safe environment.


*Ideal for a dentist to use in a mini-lab or backroom.*

### Features and Benefits:

- Small Size (11" wide x 12" deep x 8 1/2" high) fits easily on countertops.
- High Powered, reliable, 35,000 RPM V35 Electric Handpiece System ideal for most cutting and grinding needs.
- High-torque handpiece provides smooth operation and powerful performance.
- Hinged top window and internal light fixture provide bright illumination for clear viewing.
- Dual, wide hand holes allow comfortable hand positioning and access to large workspaces.
- Efficient suction fan and filter collect damaging grinding debris and maintain quiet operation.
- High capacity internal filter that is fully replaceable.
- Additional 120 V AC power outlet on back of the unit for operating electric devices.
- Hinged rear door for easy cleaning.
- Convenient handpiece clip on top of unit to keep your handpiece at the ready.



The Buffalo Dust Containment Unit  
which comes with a V35 Electric Handpiece System



### Ordering Information:

#### **Buffalo #**    **Description**

36800      *Dust Containment Unit (Illuminated Suction Cabinet and V35 Handpiece System), 120 V AC*

#### **Buffalo Dental Manufacturing Co., Inc.**

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203  
Fax: 516.496.7751 | Email: sales@mailbdm.com | www.buffalodental.com