

The Dynamic Duo^{*Patented}

Two of our most popular units combined into one complete system!



Want to save space? Save time? AND Save money?

The combination of the high-powered, reliable, V35 Electric Handpiece system and the effective, time-saving ThermaKnife will keep your lab bench organized, and keep you productive. Not only that, it will keep money in your pocket...

V35 Unit- \$563 | ThermaKnife Unit- \$245

The Dynamic Duo- \$708

That's OVER \$100.00 savings in retail pricing!



Features and Benefits:

- Four Multi-Speed options (14k, 22k, 27k, 35k rpm) allows for ideal speed selection for most cutting needs.
- High-Torque motor provides smooth operation and powerful performance with less operator fatigue.
- Fully adjustable temperature control for precise cutting action. Cutting tip heats ultra-fast to save time and eliminate waiting.
- A & B selection switch allows for the use of two different Thermal Pens to reduce time and increase trimming productivity.
- Lightweight pen features an ergonomic design and sharp, permanent tip for improved comfort and control.

Ordering Information:

Item#	Description
38760	Dynamic Duo Complete System
37590-14	Replacement Metal Chuck for Gray Handpiece
37590-37	Replacement Carbon Brushes for Gray Handpiece (2)
80510	Replacement Pen w/ Concave Tip
80520	Replacement Pen w/ Knife Edge Tip
80530	Replacement Pen w/ Waxing Tip

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203

Fax: 516.496.7751 | Email: sales@mailbdm.com | www.buffalodental.com