

# Buffalo Dental Gold Cap Burs

Quality Titanium Nitride coated carbide burs made in the USA



***Premium carbide burs provide outstanding cutting speed, durability and value!***

Using only the highest-quality materials to manufacture these burs, as well as an automated process, ensures that each and every bur produced offers the same great cutting speed, durability, and precision as the next! The premium TNT coating helps keep edges sharp and provides lubrication for superior cutting!



## Features and Benefits:

- For use on plastics, ceramics, porcelain, plaster, acrylic, stainless steel, etc. Lab burs are **NOT** for Intraoral use!
- Ultra-fast cutting and long lasting edges as a result from using only high quality materials and an automated manufacturing process
- USA manufacturing location provides and assures premium quality, while cutting out the middle men
- TNT Titanium Nitride coating provides long lasting edges and superior cutting
- Dual cut and Diamond cut heads are available. Dual cut heads improve performance and minimize clogging while Diamond cut heads help pulverize materials and excel at acrylic work

## Ordering Information:

All popular shapes, sizes and cutting styles are available. Simply contact your dealer and ask about Buffalo Dental Gold Cap Burs. If you would like more information, please contact Jeannie in Buffalo Dental Manufacturing customer service using the information below.

**Buffalo Dental Manufacturing Co., Inc.**

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203  
Fax: 516.496.7751 | Email: [sales@mailbdm.com](mailto:sales@mailbdm.com) | [www.buffalodental.com](http://www.buffalodental.com)