

M35 Compact Electric Lab Handpiece System



New High Torque Variable Speed System Features a Compact and Versatile Design

- **High Torque Cutting Power from 0 – 35,000 RPM**

New high torque handpiece and control system provide excellent cutting power and smooth control in a compact system

- **Improved Versatility**

System is designed for 120V AC operation but can be easily switched to 220V AC via setting on the underside of the Console

- **Improved Console and Controls**

Improved controls provide fully variable speed via console or foot control, simplify operation and protect the controls and handpiece

- **Excellent High and Ultra-Low Speed Operation**

System excels at providing a full range of usable operating speeds from ultra-low to 35,000 RPM

- **Heavy Duty Construction**

Components and system designed to provide reliable and trouble-free operation. One year warranty

- **Sizes**

Handpiece: 6.25 in. (159 mm) Long x 1.1 in. (28.5 mm) Back Diameter, .79 in. (20 mm) Front Diameter x 8 oz. (.23 kg) Weight

Console: 5 in. (128 mm) Width x 6.5 in. (166 mm) Depth x 3.25 in. (82 mm) Height x 2.9 lb. (1.3 kg) Weight



The Motor and Armature feature precision dynamic balancing for improved performance and less noise and vibration

- **Ordering Information**

Item #	Description
38120	M35 Electric Lab Handpiece System Complete, 120V AC (& 220V AC), including Control Console, Electric Handpiece w/HP Collet, Foot Control, Handpiece Holder and Collet Wrenches
38140	M35 Electric Lab Handpiece Only w/HP Collet Chuck
38025BS	Square Carbon Brushes for M35 Electric Handpiece, Set of 2
37805-107-3	Replacement HP Collet Chuck for M35 Handpiece, 3/32" Size
37805-COLLET3MM	Optional 3mm Collet for M35 Handpiece
37805-COLLET18	Optional 1/8" Collet for M35 Handpiece

Notes: For laboratory use only. Please use extreme caution and wear personal protective eyewear and equipment when operating the M35 system.

Buffalo Dental Manufacturing Co., Inc.,

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516-496-7200 or 800-828-0203

Fax: 516-496-7751 | Email: Sales@mailbdm.com | www.buffalodental.com

S38120 0418