

## ***Classic High Performance Vacuum Forming System in Steady and Successful Service in Thousands of Dental Offices and Labs***

### • High Performance System

Dual-stage electric vacuum motor provides powerful and immediate suction... and highly accurate adaptations- no compressed air connections required!

### • Accurate Frame and Alignment System

Spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table

### • Heavy Duty Construction

Solid construction (7 ½" W x 10 ½" D x 15" H, 32 lbs.) for maximum durability, long life and virtually no maintenance!

### • Convenient and Safe Design

Heater hood rotates behind the center post for safety and convenience, and the compound frame closure latch helps keep the components cool to the touch during repeated use


### • Wide Range of Available Accessories

Please see the complete line of Buffalo Dental vacuum forming equipment and supplies including a wide range of Vacuum Forming Sheets, Block-out Compound, Silicone Spray, Tray Handles, Knives and Burs, the ThermaKnife™ Thermal Cutting Knife, and PolyShok™ Sports Mouthguard Materials

**Made in USA!**



### • Ordering Information

<b>Item #</b>	<b>Description</b>
80186	Econo-Vac Vacuum Forming System Complete, 120V AC, incl. a reference assortment of vacuum forming materials
80187	Econo-Vac Vacuum Forming System Complete, 220V AC,  incl. a reference assortment of vacuum forming materials

*For laboratory use only.*

### **Buffalo Dental Manufacturing Co., Inc.,**

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | **Tel:** 516-496-7200 or 800-828-0203

**Fax:** 516-496-7751 | **Email:** [Sales@mailbdm.com](mailto:Sales@mailbdm.com) | [www.buffalodental.com](http://www.buffalodental.com)

S80186 0418