

Buffalo Dental Waxing Unit



Compact and Economical 3-Compartment Wax Heating Unit with Adjustable Temperature Control

Features/Benefits

- On/Off switch located on power cord to prevent wax seepage into switch mechanism
- Controlled temperatures from 140° to 210° F for convenient operation (adjustment screw on back left)
- Provides three 11/16" deep wax compartments (one is 2-1/2" x 2-7/8", and two are 1-5/8" x 1-3/8")
- Attractive design with finished metal sides, handles and helpful heat indicator light on front panel



Ordering Information

<i>Item No.</i>	<i>Description</i>
86000	Waxing Unit, 120V AC
86010	Waxing Unit, 220V AC

Buffalo Dental Manufacturing Co., Inc.,

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516-496-7200 or 800-828-0203

Fax: 516-496-7751 | Email: Sales@mailbdm.com | www.buffalodental.com

S86000 0418