

# UV-MT Model Trimmer



## Ultra-Violet Model Trimmer Kills Germs While It Turns!

*Dental stone models contain various bacteria and organisms which can be spread throughout the lab.*

*Made in USA!*

*Powerful Internal UV-LED Diode kills Bacteria and other organisms on wheel, in grindings and on model. Reduces cross contamination from model to model.*

- **Powerful and Fast Trimming**

Industry-leading extra-high torque motors and aggressive-cutting coarse Carborundum wheels cut fast and smooth saving time and money

- **Solid Construction**

Classic, heavy construction and time-proven balanced direct-drive systems run quietly and can stay steady and trouble-free on the job for decades

- **Enlarged Window**

Larger work window provides increased working area, improved vision, and full wheel use while resisting splashing

- **Versatile Table**

Calibrated work table indicates model angles and can be easily tilted and repositioned to suit nearly every trimming requirement

- **Complete Ready-to-Run System**

All systems ship complete with a pre-mounted balanced Aggressive-cutting Carborundum Trimming Wheel, ¼" ID Water Supply Tubing, 1" ID Water Exhaust Hose, and Wheel Wrenches

- **Optional Accessories**

A wide range of optional accessories are available including pre-mounted water solenoid valves, rubber feet, safety shields, numerous diamond and conventional abrasive wheels, disposable disks, and several disposable plaster trap options

- **Ordering Information**

<u>Item Number</u>	<u>Description</u>
--------------------	--------------------

61780	10" Diameter 1/3 HP Ultra-Violet Model Trimmer, 120 V AC, 60 Hz
61785	12" Diameter 1/3 HP Ultra-Violet Model Trimmer, 120 V AC, 60 Hz



**\*US and Foreign Patents Pending.**

### Buffalo Dental Manufacturing Co., Inc.,

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516-496-7200 or 800-828-0203  
Fax: 516-496-7751 | Email: [Sales@mailbdm.com](mailto:Sales@mailbdm.com) | [www.buffalodental.com](http://www.buffalodental.com)

S61780 0418