

Trap-Eze™ 2 Gallon Disposable Trap*Patented

Our new 2 Gallon Disposable Plaster Trap



It's an easy choice for under-sink waste management...



64 Oz

COMPLICATED



256 Oz

NOT SO COMPLICATED

The new Trap-Eze™ 2 Gallon Disposable Trap has been designed to fit in tighter spaces as compared to other traps on the market. Its smaller size also means a reduced overall weight when full, allowing for easy replacement by anyone!

Features and Benefits:

- Small footprint (9 1/2" H x 10" W) makes it easy to fit in smaller under-sink spaces
- Ships as a 2 pack (Regular Trap-Eze™ and refill)
- Lightweight design makes removing a full Trap-Eze™ simple and easy for anyone
- Economical trap systems designed to slow water discharge flow rate and allow sediment to settle out of waste water
- Fast and easy installation and trap replacement
- No permanent traps to clean or messy liners
- Made in the USA



Trap-Eze™ 2 Gallon,
Easy to replace!



Trap-Eze™ 2 Gallon,
Ships as 2 pack.



It would take 4 of the other brand to equal
one Trap-Eze™ 2 Gallon Trap

Ordering Information:

Item #	Description
62170	Trap-Eze™ 2 Gallon Complete Kit 2 Pack (Regular Trap-Eze™ & Refill)
62180	Trap-Eze™ 2 Gallon Refill (sealed Trap-Eze™ only), 2 pack

Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203
Fax: 516.496.7751 | Email: sales@mailbdm.com | www.buffalodental.com