

# Abbott-Robinson®

## Bristle Brushes

High quality natural bristle brushes



### **Abbott-Robinson® Bristle Brushes offer high quality and performance!**

#### **Abbott-Robinson® Bristle Brushes-**

Classic, high quality natural bristle brushes ideal for cleaning and polishing of denture acrylic, cleaning inter proximal spaces after denture de-flasking, and for general polishing applications. These brushes excel at applying polishing materials and at polishing of lab processed composites, metals, and ceramics.



No. 9 (9/16") No. 11 (11/16") No. 12 (3/4")  
Abbott-Robinson® Bristle Brushes

#### **Features and Benefits:**

- Durable design features permanently attached bristles for improved performance and extended life without flaring
- Available in Soft, Standard Stiff, and Extra Stiff bristles for use in all polishing, finishing, and cleaning tasks
- Solid Carbon-Steel plated HP mandrels for excellent fit and durability

#### **Ordering Information:**

<b>Item#</b>	<b>Description</b>		
06042	#9 X-Stiff (12 Pkgs. of 12)	06044	#9 X-Stiff (Bulk Pkg. of 144)
06012	#9 Standard-Stiff (12 Pkgs. of 12)	06014	#9 Standard-Stiff (Bulk Pkg. of 144)
05992	#9 Soft (12 Pkgs. of 12)	05994	#9 Soft (Bulk Pkg. of 144)
06102	#11 X-Stiff (12 Pkgs. of 12)	06104	#11 X-Stiff (Bulk Pkg. of 144)
06082	#11 Standard-Stiff (12 Pkgs. of 12)	06084	#11 Standard-Stiff (Bulk Pkg. of 144)
06062	#11 Soft (12 Pkgs. of 12)	06064	#11 Soft (Bulk Pkg. of 144)
06192	#12 X-Stiff (12 Pkgs. of 12)	06194	#12 X-Stiff (Bulk Pkg. of 144)
06162	#12 Standard-Stiff (12 Pkgs. of 12)	06164	#12 Standard-Stiff (Bulk Pkg. of 144)
06142	#12 Soft (12 Pkgs. of 12)	06144	#12 Soft (Bulk Pkg. of 144)

**Buffalo Dental Manufacturing Co., Inc.**

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203

Fax: 516.496.7751 | Email: sales@mailbdm.com | www.buffalodental.com